

**DEVELOPMENT OF GENOME  
EDITING CONSTRUCT TARGETING  
PHYTIC ACID AND RAFFINOSE  
FAMILY OLIGOSACCHARIDES IN  
*GLYCINE MAX*...**





A close-up photograph of several green, fuzzy soybean pods hanging from a plant. The pods are covered in fine, light-colored hairs. The background is filled with green leaves, some of which are in focus and others are blurred, creating a sense of depth. The lighting is natural, highlighting the texture of the pods and leaves.

## **INTRODUCTION:-**

- **Soybean (*Glycine max*) is a versatile and globally significant legume valued for its high protein content and numerous applications in food, feed, and industrial products.**
- **Despite its nutritional benefits, soybean contains certain antinutritional factors, such as phytic acid and raffinose family oligosaccharides (RFOs).**
- **Phytic acid can hinder mineral absorption, while RFOs can cause digestive discomfort.**
- **Understanding and regulating the genes involved in the production of these antinutritional factor is crucial for improving the nutritional quality and digestibility of soybean-based products.**

# OBJECTIVES


- **To familiarize about gene responsible for Phytic acid and Raffinose family oligosaccharide in Soybean(*Glycine Max*).**
- **To design SgRNA targeting GmIPK1a, GmIPK1b, GmIPK1c and GmGOLS1A, GmGOLS1B.**
- **Development of construct targeting genes responsible for phytic acid and raffinose family oligosaccharide.**

# IMPORTANCE AND NEED

- **Enhancing Nutritional Value:** Reducing phytic acid and raffinose family oligosaccharides (RFOs) improves mineral absorption and protein digestibility, making soybeans more nutritious ([Lin et al., 2024](#))
- **Improving Digestibility:** Lowering RFOs reduces digestive discomfort, making soy products more acceptable and easier to digest. RFOs can cause gastrointestinal issues due to their fermentation by gut bacteria([James W Anderson et al., 2009](#))
- **Supporting Human Health:** Reducing phytic acid enhances mineral uptake, preventing deficiencies and promoting better health([Anupam Roy et al., 2023](#))
- **Facilitating Food Processing:** Fewer antinutritional factors simplify processing, reduce costs, and improve product quality.
- **Advancing Genetic Research:** Understanding these genes aids in crop science advancements and can be applied to other plants for broader agricultural benefits.



# REVIEW OF LITERATURE



- In soybean seeds, PA can count for 4–5% of the seed dry weight, and 50% of the seed total phosphorus content ([Oatway et al. 2001](#)). Having six phosphate groups, PA belongs to a class of metal-chelating compounds that may bind to iron, zinc, calcium, magnesium, and potassium, among other metal cations ([Raboy 2009](#)).
- During seed germination, the PA complex releases metal cations and provides energy to support seedling growth. However, PA is considered a non-negligible anti-nutritional factor. PA lowers the absorption and utilization of minerals (e.g., iron, zinc, calcium) and negatively affects human health for people subsisting on grain-based diets, mainly in developing countries ([Zhou and Erdman 1995](#)).
- The major genes responsible for the low phytic acid germplasm include the inositol pentose-phosphate kinase gene (GmIPK1a) in PA biosynthesis ([Yuan et al. 2012](#)), and the vital multi-drug-resistant protein 5 gene (GmMRP5a) involved in PA transport ([Gillman et al. 2009](#)).


# REVIEW OF LITERATURE

- Studies explored the use of CRISPR/Cas9 technology for gene editing in soybean, specifically targeting two genes involved in the phytic acid synthesis pathway: GmIPK1 and GmIPK2([Carrijo et al., 2021](#)).
- GmIPK1 and GmIPK2 were chosen as target genes due to their roles in phytic acid biosynthesis in soybean seeds. GmIPK1 is involved in inositol-pentakisphosphate 2-kinase activity, while GmIPK2 likely contributes to the biosynthesis of phytic acid ([Carrijo et al., 2021](#)).
- The successful application of CRISPR/Cas9 technology in editing genes related to phytic acid biosynthesis in soybean highlights the potential for targeted crop improvement and trait modification through genetic engineering ([Carrijo et al., 2021](#)).

# REVIEW OF LITERATURE

- Results in studies showed that higher-order mutant lines with multiplex mutations had significantly decreased seed PA content. However, these mutants exhibited poor agronomic performance([Lin et al., 2023](#)).
- In contrast, two lines with single mutations in IPK1b and IPK1c showed moderately reduced PA content and regular agronomic performance compared to the wild type ([Lin et al., 2023](#)).
- This suggests that moderately decreasing PA by targeting single GmIPK1 genes, rather than multiplex mutagenesis for ultra-low PA, is an optimal strategy for low-PA soybean with minimal trade-off in yield performance ([Lin et al., 2023](#)).

# REVIEW OF LITERATURE



- The raffinose family oligosaccharides (RFOs) are a major class of water-soluble carbohydrates present in soybean seeds. Soybean RFOs include raffinose, stachyose and verbascose which, respectively, are mono-, di-, and tri-galactosides of sucrose ([Kennedy et al., 1985](#); [Hou et al., 2009](#)), with stachyose being the most prominent RFO in mature seeds.
- GOLS catalyzes the first step in the RFO biosynthetic pathway ([Nishizawa et al., 2008](#)). Therefore, study hypothesized that knock-out mutations in the *GmGOLS* genes, specifically in *GmGOLS1A* that encodes the major GOLS expressed in developing soybean seeds, would affect the final sugar composition in mature seeds.
- RFOs were proposed to be ubiquitous storage products and factors of desiccation tolerance in soybean seeds, with a potential role in other stress responses ([Taji et al., 2002](#); [Wang et al., 2009](#); [Loedolff, 2015](#); [de Souza Vidigal et al., 2016](#)). However, previous studies of soybean seed lines with low stachyose and/or raffinose showed no significant differences in stress tolerance or germination compared to wild-type ([Neus et al., 2005](#); [Valentine et al., 2017](#)). RFOs are indigestible by human and other monogastric animals that do not produce  $\alpha$ -galactosidase to break the  $\alpha$ -1,6-glycosidic bonds in RFOs.



# REVIEW OF LITERATURE

- Studies aimed to investigate the impact of knocking out the GmGOLS1A (Glyma.03G222000) and GmGOLS1B (Glyma.19G219100) genes, which encode Galactinol Synthase, on Raffinose Family Oligosaccharide biosynthesis in soybean seeds([Le et al., 2020](#)).
- The researchers employed a CRISPR/Cas9 system with dual guide RNAs to induce targeted knockouts in the GmGOLS1A and GmGOLS1B genes. Genotyping of T0 plants confirmed efficient disruption of the target sites.([Le et al., 2020](#))
- The knockout of GmGOLS1A and GmGOLS1B genes led to a significant decrease in the total Raffinose Family Oligosaccharide content in soybean seeds compared to wild-type seeds.([Le et al., 2020](#))

# MATERIALS AND METHODS

## *In Silico* work:-

- **Designing of desired sgRNA with the use of sgRNA designing tools such as CRISPR-P & CRISPR-GE.**
- **Sequence of plasmid pDIRECT22\_A (from Addgene).**
- **Using sequence in SNAP GENE tool.**
- **Annealing Oligos and Cloning by inserting fragments of vector(AaR1) and Annealed Oligose.**
- ***In silico* PCR for cloned DNA fragment with Primers ( gRNA and NB463 sequences)**



## ***In vitro* work:-**

- **Formation of sgRNA oligoduplex.**
- **Materials which we will need for this reaction will be:-**


<b>Chemicals</b>	<b>Quantity</b>
<b>Forward Primer(100μM)</b>	<b>3μl</b>
<b>Reverse Primer(100μM)</b>	<b>3μl</b>
<b>T4 DNA Ligase Buffer</b>	<b>3μl</b>
<b>T4 Polynucleotide kinase</b>	<b>2μl</b>
<b>MQ water</b>	<b>19μl</b>

- **Followed by PCR:-**  
**37°C (60 min)**  
**95 °C (5 min)**  
**25 °C (0.1 °C/sec RAMP mode or water bath)**
- **Isolation of Plasmid(pDIRECT22\_A).**

## **GOLDEN GATE REACTION:-**

- **Involves simultaneously Restriction digestion and Ligation.**
- **We will need:-**

Chemicals	Quantity
pDIRECT_22A	50ng/μl
25x diluted oligo duplex	1μl
Restriction enzyme AarI	0.5μl
AarI Oligonucleotide	0.4-0.5μl
T4 DNA Ligase Buffer	2μl
T4 DNA Ligase	1μl
MQ water	depends
Final volume	20μl



**Followed by PCR:- 37°C(5 min), 16°C(10min), 37°C(15 min), 80°C(5 min, for inactivation)**



# **PREPARATION OF ULTRACOMPETENT CELLS**

## **Materials and Reagents:-**

- 1. E. coli strain (e.g., DH5 $\alpha$ )**
- 2.LB medium (Luria-Bertani)**
- 3.LB plates with the appropriate antibiotic**
- 4.Transformation buffer:**
  - 1. PIPES (1 M, pH 6.7)**
  - 2. CaCl<sub>2</sub> (0.1 M)**
  - 3. KCl (1.5 M)**
  - 4. MnCl<sub>2</sub> (0.55 M)**
- 5.Sterile distilled water**
- 6.Glycerol (10%)**
- 7.Sterile microcentrifuge tubes**
- 8.Ice-cold 100 mM CaCl<sub>2</sub> solution**

**PROCEDURE OF MAKING COMPETENT CELLS of *E. coli* :-**

**•Preparation of Inoue Transformation Buffer:-**

- 1. Dissolve 30.24 g of PIPES in about 80 mL of sterile distilled water.**
- 2. Adjust the pH to 6.7 with KOH.**
- 3. Add 10.9 g of  $\text{CaCl}_2 \cdot 2\text{H}_2\text{O}$ , 37.2 g of KCl, and 4.95 g of  $\text{MnCl}_2 \cdot 4\text{H}_2\text{O}$ .**
- 4. Adjust the final volume to 100 mL with sterile distilled water.**
- 5. Sterilize the solution by passing it through a 0.22  $\mu\text{m}$  filter.**
- 6. Store the buffer at 4°C. Do not autoclave.**

**• Growth of Bacterial Culture:-**

- 1. Inoculate a single colony of *E. coli* into 50 mL of LB medium.**
- 2. Incubate the culture overnight at 18-22°C with shaking at 250-300 rpm.**
- 3. The next day, dilute the overnight culture 1:100 into fresh LB medium.**
- 4. Incubate at 18-22°C with shaking until the  $\text{OD}_{600}$  reaches 0.55-0.6 (mid-log phase).**



**•Preparation of Competent Cells:-**

- 1. Chill the culture on ice for 10 minutes.**
- 2. Transfer the culture into pre-chilled sterile centrifuge bottles.**
- 3. Centrifuge at 2500 x g for 10 minutes at 4°C.**
- 4. Discard the supernatant and gently resuspend the pellet in 80 mL of ice-cold Inoue transformation buffer.**
- 5. Incubate on ice for 10 minutes.**
- 6. Centrifuge again at 2500 x g for 10 minutes at 4°C.**
- 7. Discard the supernatant and gently resuspend the pellet in 20 mL of ice-cold Inoue transformation buffer.**

**•Addition of Glycerol:-**

- 1. Add 1.5 mL of sterile 100% glycerol to the cell suspension to a final concentration of 7%.**
- 2. Mix gently by swirling.**

**•Aliquoting and Storage:-**

- 1. Aliquot 200 µL of the competent cell suspension into pre-chilled sterile microcentrifuge tubes.**
- 2. Flash freeze the aliquots in liquid nitrogen.**
- 3. Store the frozen aliquots at -80°C.**

***E. Coli* Transformation procedure:-**

- 1. Keep the prepared competent cells in ice until Thaw.**
- 2. Add 250ng of complete GOLDEN GATE REACTION product.**
- 3. Tap Gently.**
- 4. Keep competent cells in ice for 15 min.**
- 5. Transfer vial 42°C in water bath/ heating block for 90 secs.**
- 6. Transfer vial to ice.**
- 7. Add 1ml of LB Broth to vial.**
- 8. Incubate the vial culture at 37°C with shaking for 1.5 hrs.**
- 9. After incubation centrifuge the culture at 5000 rpm for 5 minutes.**
- 10. Discard the supernatant (precaution; let some LB broth be inside).**
- 11. Resuspend the cell pellet with LB broth (100µl)**
- 12. Spread the culture on selection plate.**
- 13. Spread the aliquot on petri plate until it become completely dry(LB plates containing X Gal, IPTG and Kanamycin).**
- 14. Incubate the plate at 37 °C for 12 hrs.**
- 15. Select white colonies (Blue white screening method) and do Colony PCR Forward Primer(sgRNA) and Reverse Primer(NB463).**
- 16. If we get 429bp amplicon as a result then we will send it for sequencing for further confirmation.**



# PLAN OF WORK

- Identification of target sequence
- sgRNA designing
- Primer designing and Oligoes designing
- *In silico* development of CRISPR construct
- Isolation of Plasmid
- DNA oligo duplex preparation
- Cloning (Golden Gate Reaction)
- Competent cell preparation(E Coli DH5α)
- *E Coli* Transformation(Heat and Shock)
- Confirmation and Validation of plasmid(PCR)
- Agrobacterium Transformation( Freeze Thaw method)

# IMPLICATIONS

- **Knowledge of CRISPR Cas9, gRNA design, competent cell preparation, cloning, Transformation.**
- **Proficiency in tools like CRISPR-P, CRISPR-GE, CHOP-CHOP, Soykb, Phytozome, Soybase.**
- **Proficiency in softwares like SNAPGENE for designing, analysis, visualization.**
- **Crop improvement, sustainable agriculture, ethical implications.**
- **Competent cell preparation, cloning and transformation, gRNA design knowledge.**
- **Gene editing knowledge will help in doing something beneficial for mankind.**

# ACTIVITIES PROPOSED AND WORK SCHEDULE

Task	Duration
Identification of desired genes and gene IDs	3 days
Valid gRNA design	3 days
Primers/Oligoduplex in silico formation	1 day
Snapgene	1 week
DNA oligoduplex formation (in vitro)	1 day
Golden Gate reaction (Digestion and Ligation)	1 day
Preparation of competent cells	2 days
E. coli Transformation	1 day
Sequencing and Validation	2 weeks
Agrobacterium Transformation	1 day
Validation of Transformation	1 day



# REFERENCES

Lin, W., Bai, M., Peng, C. *et al.* Genome editing toward biofortified soybean with minimal trade-off between low phytic acid and yield. *aBIOTECH* (2024). <https://doi.org/10.1007/s42994-024-00158-4>

Le, Huy, et al. "CRISPR/Cas9-mediated knockout of galactinol synthase-encoding genes reduces raffinose family oligosaccharide levels in soybean seeds." *Frontiers in Plant Science* 11 (2020): [612942](https://doi.org/10.3389/fpls.2020.612942).

Carrijo, J., Illa-Berenguer, E., LaFayette, P., Torres, N., Aragão, F.J., Parrott, W. and Vianna, G.R., 2021. Two efficient CRISPR/Cas9 systems for gene editing in soybean. *Transgenic Research*, 30(3), pp.239-249.

James W Anderson, Pat Baird, Richard H Davis, Stefanie Ferreri, Mary Knudtson, Ashraf Koraym, Valerie Waters, Christine L Williams, Health benefits of dietary fiber, *Nutrition Reviews*, Volume 67, Issue 4, 1 April 2009, Pages 188–205, <https://doi.org/10.1111/j.1753-4887.2009.00189.x>

Kumari A, Roy A. Enhancing micronutrient absorption through simultaneous fortification and phytic acid degradation. Food Sci Biotechnol. 2023 Jan 27;32(9):1235-1256. [doi: 10.1007/s10068-023-01255-8](https://doi.org/10.1007/s10068-023-01255-8). PMID: 37362807; PMCID: PMC10290024.

*Thank  
You*